




## Preat O-Ring Size Chart



To determine the correct size O-ring needed, measure the retention zone (red line) – or the area on the male post where the O-ring seats. Select the O-ring from the size chart with the closest matching internal diameter to assure proper fit. The O-ring Combo pack is a great option as well.



O-RING REPLACEMENT SIZE CHART				
		Internal Diameter (ID)	Outside Diameter (OD)	
	SIZING	mm	mm	PRICE
Black #2 Pkg 12	Micro	1.00 mm	3.60 mm	\$27.00
Black #2M Pkg 6	1.8 Mini	1.20 mm	3.60 mm	\$27.00
White #3 Pkg 12	Standard	1.40 mm	4.50 mm	\$23.00
Black #3 Pkg 12	Standard	1.40 mm	4.50 mm	\$27.00
Black #4 Pkg 6	Large	2.00 mm	5.20 mm	\$27.00
Black #5 Pkg 6	Extra Large	2.40 mm	6.25 mm	\$30.00
Combo Pack 2 of each ring (10 total)	Micro (2) Mini Implant (2) Standard (2) Large (2) Extra large (2)	2 of each size (Black)		\$48.00
Black #2 Pkg 12	Micro (2) Mini Implant (2) Standard (2) Large (2) Extra large (2) Insertion Tool	2 of each size (Black) Insertion Tool		\$100.00
		Metal Retainer Ring #2 (fits both #2 and #2M) Metal Retainer Ring #3 (each)		\$18.00
		Metal Retainer Rings #4 (each)		\$27.00
		Metal Retainer Rings #5 (each)		\$27.00

Contact Preat at 1-800-232-7732 or visit [preat.com](http://preat.com) to order online